



Developability Assessment Services

LakePharma Hayward

NR-7041.20191223

About LakePharma

- LakePharma is a US-based Biologics CRDMO with operations in CA, MA, and TX
- Specializes in the production and evaluation of DNA vectors, viral vectors, cell lines, proteins, antibodies and conjugates, while providing integrated solutions bridging discovery, engineering, development, and GMP manufacturing. *Ask us about ISPs!*
- LakePharma has contributed to the development of 200+ therapeutic or diagnostic products and strives to do hundreds more.
- LakePharma Hayward is ISO 9001:2015 certified, and LakePharma Hopkinton is ISO 13485:2016 certified.



LakePharma Hayward


Dedicated facility for upstream and downstream process development, research cell bank (RCB) generation, and GMP manufacturing of analyte specific reagents.

- All developability assessments are performed at this site



Developability: Assess the Path to Product

LakePharma offers developability studies for antibody drug candidates with promising biochemical and biophysical properties.



LakePharma's Developability Assessment Packages

Rapid and small-scale assessment of drug candidates

Developability Package 1

- *In silico* Sequence Liability analysis
- *In silico* Immunogenicity analysis
- Turnaround time: 1 week

Developability Package 2

- Polyspecificity Assessment
- Integrity and Stability Assessment
 - Aggregation
 - Purity
 - Charge Variant
 - Thermostability
 - Post-translational modifications
- Turnaround time: 2-3 weeks

Formulation and stability study

Developability Package 3

- Buffer Exchange
 - LakePharma standard panel formulations
 - Client may opt to choose their buffers
- Forced Degradation
 - Thermal stress
 - Freeze thaw
- Integrity and Stability Assessment following stress conditions
- Additional available stress conditions
 - Agitation
 - Oxidation
 - Photostability (Light)
 - pH acid/base

Therapeutic Developability Assessment Packages 1 and 2



Developability Package 1

In silico predictive tools are employed to serve as selection criteria for better and safer therapeutic leads.

1. *In silico* sequence liability analysis
2. *In silico* immunogenicity analysis



Developability Package 2

A series of fast and small-scale bioanalytical tests will be performed to determine the feasibility and liability of the drug candidates prior to drug development.

1. Polyspecificity assessment
2. Integrity and stability assessment

Key Features:

- Rapid and small-scale assessments for drug candidates
 - *In silico* analysis can be completed in 1 week
 - Polyspecificity and integrity assessment can be completed in 2-3 weeks
- Only a small amount of materials is needed for stability and liability determination

Therapeutic Developability Assessment Package 3

The goal of this package is to assess the stability of the molecule and identify a suitable formulation for the antibody.



Developability Package 3

1. Buffer Exchange

- LakePharma standard panel formulation
- Client may opt to choose their buffers

2. Forced Degradation

- Thermal stress
- Freeze thaw

3. Integrity and Stability Assessment

- Aggregation
- Purity
- Charge Variant
- Thermostability
- Post-translational modifications

4. Additional available stress conditions

- Agitation
- Oxidation
- Photostability (Light)
- pH acid/base

Key Features:

- Rapid formulation studies for drug candidates
 - 2 weeks incubation time for standard forced degradation
 - Formulation studies can be completed in 5 - 6 weeks
- Robust formulation assessments with minimal materials
 - LakePharma standard panel includes common formulations used for commercial antibodies
 - Various bioanalytical and stability tests are provided



Developability Assessment Checklist

- ☐ *In silico* sequence liability analysis and immunogenicity analysis
 - Heavy and light chain variable liability analysis
 - VH and VL analysis
 - Nine core residues with potential affinity to MHC class II identified




- ☐ Productivity readiness check
 - Transient production in HEK293 or CHO system

- ☐ Integrity and stability check (to demonstrate whether the antibody is stable biochemically)
 - Intact mass/peptide mapping/ PTM by mass spec
 - DSF/DLS – thermostability assessment
 - Aggregation/fragmentation/PTM and glycan profiling over 2 weeks incubation under stressed conditions
 - SE-HPLC/LabChip CE-SDS


- ☐ PK readiness check
 - Poly-specificity ELISA, surface hydrophobicity assay – specificity requirement




DSF/SLS



HPLC



Poly-Specificity ELISA




CASE STUDY

Assessing a mAb Candidate Utilizing Package 2



A Therapeutic mAb Candidate Ready to Move Forward



Purity requirement **rCE>95%**

The mAb is ready to move forward as results show >95% purity by reduced CE-SDS analysis and >90% purity by non-reduced CE-SDS analysis

Purity requirement **nrCE>90%**


Sample ID	%LMWC	%(HC+LC)	%MMWC	%HMWC
mAb reduced	N/A	99.3	0.7	N/A
mAb non-reduced	6.5	93.5	N/A	N/A



LabChip® GXII Touch HT

A Therapeutic mAb Candidate Ready to Move Forward

SE-HPLC – monomeric state



The mAb is ready to move forward as results show >90% peak area for the monomer peak.


Peak #	Peak Size (kDa)	Peak Area %	Peak ID
1	~1500	1.4	Aggregate
2	~250	2.3	Aggregate
3	~150	96.3	Monomer



Waters® ACQUITY UPLC® System

A Therapeutic mAb Candidate Ready to Move Forward

Intact mass measurement – integrity of the molecule



The mAb is ready to move forward based on intact mass data. The heavy and light chain mass was confirmed to match the sequence's expected mass differences.


	Heavy Chain	Light Chain
Measured Mass	48636.90 Da	23794.40 Da
Calculated Mass	48634.90 Da	23794.53 Da
Delta Mass	+2.00 Da	-0.13 Da



Xevo™ G2-XS QToF
Mass Spec System

A Therapeutic mAb Candidate Ready to Move Forward

cIEF – Charge variants profiling



Peak Notes	pI	Peak Area %	# Peak ID
1-2	---	15.38	Acidic Peak Group
3	8.28	65.28	Main Peak
4-5	---	19.34	Basic Peak Group


The mAb is ready to move forward as two pI markers are observed and a charge profile of the mAb is produced with a main peak observed



iCE3™ system

A Therapeutic mAb Candidate Ready to Move Forward

DSC – Thermostability





Peak #	Tm (°C)	ΔHcal (kJ/mol)
1	71.3	438
2	80.7	1487
3	82.1	808



Nano DSC

The mAb is ready to move forward as melting temperatures of 71.3°C, 80.7°C, and 82.1°C were observed. The three peaks signify the Fc (C_{H2}, C_{H3}) and Fab regions of the mAb.


DSF – Thermostability and Aggregation onset



UNcle

A Therapeutic mAbCandidate Ready to Move Forward

Affinity measurement against antigen in BLI or SPR – affinity requirement



The mAb is ready to move forward as a binding affinity of $1.3\text{E-}11$ was observed


Method	Loading Sample ID	Sample ID	KD (M)	kon(1/Ms)	kdis(1/s)	Full X ²	Full R ²
Octet	Antigen	mAb #1	$1.3\text{E-}11$	$1.6\text{E+}05$	$2.0\text{E-}06$	0.0217	0.9998



Octet®

Commercial mAb Polyspecificity Assessment

Polyspecificity ELISA– Specificity requirement



HEK293 whole cell lysate (WCL) and BVP ELISA showed consistent results, which serve as good PK indicators.

Working with LakePharma

- Comprehensive technology platform
- Technical consultation with experts in antibody discovery, protein production, and GMP manufacturing
- Online data system for 24-hour access to project information (timelines, data, team communications)
- Strong project management with regular project updates (email and teleconferences)

For more information, please contact Inquiries@LakePharma.com